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
pUC19 DNA

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
Instructions for use

pUC19 DNA

High-purity plasmid DNA from pUC19 (pUC19c, Acc. No.: L09137) with a total length of 2686 base pairs.

- oc-Form
- ccc-Form

1,2 % Agarosegel

OD_{260/280}: >1,70
genomic DNA: <2 %
oc-form: <3 %
ccc-form: >95 %
base pairs: 2686


Design: (gene name)

Lac-Operon 682-237: (*lacZ*α) :*lacZ*
Ampicillin resistance (*bla*): *amp* / 2417-1629
Origin of replication (*ori*): *ori* / 1455-870

The coordinates are listed clockwise in the plasmid map.

Description

pUC19 is a standard high-copy cloning vector for *E.coli* recombinants. The position of the Multi Cloning Site (MCS), which is inserted in frame into the *LacZ*-gene, enables a blue-white selection of insert containing plasmid DNA by α-complementation. The copy number of plasmids per cell depends on the temperature and equals app. 70-80 at 37 °C and over 200 at 42 °C. The total sequence can be obtained from the EMBL data bank (ebi.ac.uk/embl/) under accession number (L09137).

Multi-Cloning-Site:

5'-ATG ACC ATG ATT ACG CCA AGC TTG CAT GCC TGC AGG TCG ACT CTA GAG GAT CCC CGG GTA CCG AGC TCG AAT TCA CTG GCC-3'

Hind III Sph I Pst I Sal I Xba I BamH I Sma I Kpn I Sst I EcoRI
AccI Hinc II Xma I

The **pUC19 plasmid-DNA** was produced in host bacteria low in nuclease, and isolated through ion-exchange chromatography. It is available at an extremely high percentage in "supercoiled" form (ccc-form >95 %) and being a high-purity DNA, (OD_{260/280} >1.70) it is ideal for different applications:

- controlling transformation efficiency
- comparison with plasmid-ccc-forms of unknown size
- cloning DNA-fragments
- restriction (production of DNA-linear markers)
- calibrating DNA-precipitating or DNA-preparing procedures or equipment
- blocking of DNA-chips with plasmid-DNA of a defined sequence.

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Carl Roth GmbH + Co. KG