ROTI®DipSlides PCA/VRBD

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Product Data Sheet

ROTI®DipSlides PCA/VRBD

DipSlides and contact agar in one versatile tool. Sterile nutrient broth for total germ count and for detection of enterobacteria
2924

- ✓ Ready-to-use and sterile
- ✓ Very easy application, reliable assay and easy interpretation
- ✓ Suitable for analysis of each surface and a multitude of liquids like water, milk, water-soluble paints, cosmetics, urine, etc.
- ✓ For microbiological laboratories, clinics and for hygiene control in food production

COMPOSITION

PCA (side A): Plate Count Agar (acc. to APHA, ISO 4833:2003) (X930) for total germ count (**amber**)

VRBD (side B): Violet Red Bile Agar with Glucose (Ph. Eur.) (X939) for detection of enterobacteria. (**purple-red**)

STORAGE, EXPIRY AND SHIPMENT

ROTI®DipSlides are shipped at room temperature. Please make sure that ROTI®DipSlides are stored at +15 – 25 °C. Store boxes lid-up, to avoid moistening of the agar with condensed water.

ROTI®DipSlides may not be frozen. Slides that have been frozen have to be disposed of, since the medium is likely to have taken severe damage and growth of microorganisms is significantly hindered.

Store at +15 – 25 °C and protected from direct light.

Expiry of ROTI®DipSlides is 5 months after production. Expiry date is noted on each package and on each vial with „EXP“.

COMPOSITION OF AGARS (in g/l, approximated)

For further information please see the respective data sheets of the agars.

PCA (X930)	Casein Peptone 5.0, Glucose 1.0, Yeast Extract 2.5, Agar 15.0, Final pH 7.2 ± 0.2. Colour: amber	VRBD (X939)	Pancreas Hydrolysate from Gelatin 7.0, Yeast Extract 3.0, Glucose-monohydrate 10.0, Sodium Chloride 5.0, Bile Salts 1.5, Neutral Red 0.03, Crystal Violet 0.002, Agar 15.0, Final pH 7.4 ± 0.2. Colour: purple-red
---------------	--	----------------	--

APPLICATION

1. Remove carrier from the sterile vial. The screw cap serves automatically as holder.
2. Press agar gently onto the surface that has to be tested. Alternatively, immerse the agar for several seconds in the liquid, in order to wet both sides completely. Carefully allow excess sample to run off. In case of samples of high viscosity, we recommend diluting the sample by 1:2 to 1:10. Dilution must be taken into account when calculating the sample contamination. Please note: In case the DipSlides are used as contact agar, two different surfaces may be tested. We strongly recommend, however, NOT to use one DipSlide for testing two different liquids, since the risk of cross-contamination is very high! Do NOT touch the agar (risk of contamination!)
3. Replace carrier into the vial and close the cap moderately.
4. Incubate vial with carrier at a temperature suitable for the bacteria (e.g. 35 ± 2 °C for 18 - 24 h). During incubation, the vials should be screwed on only loosely, in order to avoid condensation of water. Please note: In accordance to official regulations or internal protocols other incubation parameters may be chosen.
5. After the appropriate incubation time span count and register colonies.

EVALUATION OF RESULTS

After the incubation time span required colonies are counted and noted as cfu (colony forming units). For this purpose, the carrier should not be removed from the vial, in order to prevent release of the microorganisms. After incubation, tightly close the tube and let it cool for at least 2 hours at 4 °C. Condensed water may then be tapped into the lid, enabling detailed analysis of the agars through the plastic wall.

Any validation required has to be undertaken by the user. We recommend to test each sample liquid or surface at least twice (with two ROTI®DipSlides), in order to validate results. While calculating sample contamination, eventual sample dilution factors have to be taken into account. Evaluation of results is done semi-quantitatively according to the table given below.

For solutions the following approximation may be applied:

Cfu of one side x 40 is approx. cfu / ml solution. Detection limit is 10² CFU / ml solution.

Table 1: Evaluation

Colony no.	cfu / ml *	Result	Interpretation
0	<10 ²	Negative	Contamination below lower detection limits
1-5	ca. 10 ²	Weakly positive	Weak contamination of solution/surface
10-50	ca. 10 ³	Positive	Low contamination of solution/surface
100-500	ca. 10 ⁴	Strongly positive	High contamination of solution/surface
>500	>10 ⁴	Highly positive	Very high contamination of solution/surface

* This calculation may be used as approximation to numbers only. For definite quantitation of contaminated solutions or surfaces, sophisticated microbiological procedures have to be performed using Roth nutrient dry media.

DISPOSAL

Disposal of ROTI®DipSlides bearing cultures has to be done according to all common regulations for disposal of the respective microorganisms, and to regulations of good laboratory practice. E.g. autoclaving for 20 mins. at 121 °C or incubation in 70 % ethanol.

MICROBIOLOGICAL TEST

PCA (35 ± 2 °C / 18-24 h bzw. *: 30 °C / 72 ± 3 h)

Test Strain	Growth	Colour
* <i>Escherichia coli</i> ATCC 25922	Good	light-beige
<i>Escherichia coli</i> ATCC 8739	Good	light-beige
<i>Staphylococcus aureus</i> ATCC 25923	Good	light-beige
* <i>Staphylococcus aureus</i> ATCC 6538	Good	light-beige
<i>Staphylococcus epidermidis</i> ATCC 12228	Good	light-beige
* <i>Bacillus subtilis</i> ATCC 6633	Good	light-beige

VRBD (30-35 °C / 18-24 h and *: 37 °C / 24 ± 2 h, respectively)

Test Strain	Growth	Colour
* <i>Escherichia coli</i> ATCC 11775	Good	Red
* <i>Escherichia coli</i> ATCC 8739	Good	Red
* <i>Salmonella gallinarum</i> NCTC 9240	Good	Red
<i>Staphylococcus aureus</i> ATCC 6538	Inhibited	-
* <i>Shigella flexneri</i> ATCC 29903	Good	Red
<i>Pseudomonas aeruginosa</i> ATCC 9027	Good	-
<i>Enterococcus faecalis</i> ATCC 29122	Inhibited	-

PLEASE NOTE!

Application and evaluation of results of ROTI®DipSlides may only be carried out by trained personnel. The current instructions-for-use gives only general information for test performances, which have to be adopted by each user for his or her special requirements, particularly if specific regulations have to be followed. Any validation required has to be undertaken by the user.

Table 2: Overview

Ord. No.	Product	Agar Side A	Agar Side B	Application	Inc. temp.	Inc. time
2924.1	ROTI®DipSlide PCA/VRBD	Plate Count (X930)	Violet Red Bile with Glucose (X939)	whole germs + enterobacteria	35 ± 2 °C	18 - 24 h
2925.1	ROTI®DipSlide PCA/VRB	Plate Count (X930)	Violet Red Bile with Lactose (X940)	whole germs + coliforms	35 ± 2 °C	18 - 24 h
2926.1	ROTI®DipSlide PCA/RBCplus	Plate Count (X930)	Bengal Red with chloroamph. (AE24)	whole germs + yeast/moulds	20 - 25 °C	48 h* / 5 d
3934.1	ROTI®DipSlide CASO/RBCenr	Tryptic Soy Agar (CASO) (X937) with TTC	Bengal Red with antibiotics	whole germs + yeast/moulds	25 - 30 °C	48 h* / 5 d
2927.1	ROTI®DipSlide PCA	Plate Count (X930)	Plate Count (X930)	whole germs	35 ± 2 °C	18 - 24 h
2954.1	ROTI®DipSlide VRBD	Violet Red Bile with Glucose (X939)	Violet Red Bile with Glucose (X939)	enterobacteria	35 ± 2 °C	18 - 24 h
2956.1	ROTI®DipSlide Sab4plus	Sabouraud 4 % Glucose with chloroamph.(X933)	Sabouraud 4 % Glucose with chloroamph.(X933)	yeast/moulds	30 ± 2 °C	48 h - 7 d
2957.1	ROTI®DipSlide CLED/MacCon	CLED (2835)	MacConkey (X922)	urinary tract bacteria + enterobacteria/ coliforms	35 ± 2 °C	18 - 24 h

* Register PCA and CASO-TTC results after 48 h.

ROTI®DipSlide PCA/VRBD

20 slides

2924.1

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 07/2021