

Specification

Item number:	8356		
1-O-Methyl-β-D- glucopyranoside			
$\geq 99\%$, for biochemistry			
CAS No.:	709-50-2	print date:	22.04.2017
formula:	$C_7H_{14}O_6$		
density:		storage temperature:	+4 °C
molecular weight:	194,19 g/mol		

Type analysis

Appearance	white to off-white powder
Assay (HPLC)	$\geq 99,0\%$
Solution (1 % in H₂O)	colourless
Specific rotation $[\alpha]_D^{20}$ (1 % in H₂O, dry substance)	-35,0° to -31,0°
Water (KF)	$\leq 6,0\%$

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

U. Jeckel
Quality Assurance

S.Hildebrandt
Quality Management

This document is computer printout and has therefore not been signed by hand.