

Specification

Item number: HP67

Agarose-Tablets

Roti@garose 0,5 g agarose/tablet

CAS No.: 9012-36-6

print date: 07.11.2016

formula:

density:

molecular weight:

Guarantee analysis

Electroendosmosis (EEO)	0,05-0,1
Water content	<10 %
Sulphate (SO ₄)	<0,2 %
Gel strength (1.5 % gel)	≥1200 g/cm ²
Gelling temperature	37 ± 2,0 °C
Melting temperature	89 ± 1,0 °C
DNases, RNases	none detected

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

O. Lautenschläger
Head of Quality Control

S.Hildebrandt
Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Vorsitzender des Aufsichtsrats: Eberhard Gaul
Geschäftsführer: Dr. Alfred Wagner und Lothar Haidmann

