

Specification

Item number:	8763		
Betulin(ol)			
ROTICHROM® ~CHR			
CAS No.:	473-98-3	print date:	21.04.2017
formula:	$C_{30}H_{50}O_2$		
density:		storage temperature:	+4 °C
molecular weight:	442,70 g/mol		

Appereance	white to beige Powder
DC	1 Main Spot, 1 Side Spot
Mobile Phase	Ethyl acetate:n-heptane(50:50; V/V)
Stationary Phase	Silica gel 60 F 254, thickness 0,2 mm
Solvent	Ethylacetate
Detector	Anisic aldehyde-sulphric acid reagent

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

U. Jeckel
Quality Assurance

S.Hildebrandt
Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Vorsitzender des Aufsichtsrats: Eberhard Gaul
Geschäftsführer: Lothar Haidmann

