

Bromičnan draselný

Potassium bromate

KBrO₃

M_r 167,01

CAS: 7758-01-2

EINECS: 231-829-8

11090

p.a. /A.G.

Obsah / Assay..... min. 99,0%
pH (5% v H₂O, 20 °C) / (5% in H₂O, 20°C)..... 5,0 - 9,0
Bromidy / Bromide max. 0,05%
Sířany / Sulfate..... max. 0,01%
Železo / Iron max. 0,005%
Olovo / Lead max. 0,005%

11080

čistý / pure

Obsah / Assay..... min. 99,0%

Klasifikace podle nařizení 1272/2008/ES / Classification according to 1272/2008/ES



H 271, 350, 301

P 201, 309+311, 210, 221

Signální slovo: Nebezpečí / Signal word: Danger

Klasifikace podle ADR/RID / Classification according to ADR/RID

ADR/RID: 5.1/ PG II
UN 1484