

Specification

Item number:	2247		
Dextran 60			
for biochemistry			
CAS No.:	9004-54-0	print date:	01.02.2019
formula:	$(C_6H_{10}O_5)_n$		
density:			
molecular weight:	~ 60000 g/mol		

Type analysis

Appearance	white to off-white powder
Molecular mass	50000-76000
pH-value	4,5-7,0
Specific rotation $[\alpha]_D^{20}$	+180° to +201°
Loss on drying	≤ 7 %

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

U. Jeckel
Quality Assurance

S.Hildebrandt
Quality Management

This document is computer printout and has therefore not been signed by hand.