

Specification

Item number:	1C33		
4-Chlorophenol			
≥98 %, for synthesis			
CAS No.:	106-48-9	print date:	08.05.2021
formula:	C ₆ H ₄ ClOH		
density:	1,4		
molecular weight:	128,56 g/mol		

Type analysis

Appearance	colourless to slightly brown crystals or fused mass
Assay (GC)	≥98 %
Melting range	39-45 °C

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

Dr. R. Niemand
Head of Quality Assurance

S. Lorsee
Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Vorsitzender des Aufsichtsrats: Eberhard Gaul
Geschäftsführer: André Houdelet

