

Specification

Item number: 3727

RNAse Inhibitor

40 U/μl, for biochemistry and molecular biology

CAS No.:

print date:

01.01.2025

formula:

density: 1,18

storage temperature: -20 °C

molecular weight:

Guarantee analysis

Concentration	40 U/μl
Purity	≥90 %
RNAse contamination	passed
Endonuklease Contamination	passed
Exonuklease Contamination	passed

This RNAse inhibitor inhibits eukaryotic type ribonucleases e.g.

RNAse A, B and C.

Procaryotic RNases such as RNAse 1, RNAse T1, RNAse T2, S1 nucleases and RNAse H are not inhibited.

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

Dr. R. Niemand

Head of Quality Assurance

S. Lorsee

Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5

76185 Karlsruhe

Telefon 0721/5606-0

Telefax 0721/5606-149

E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

