

Specification

Item number: 1A9C

Acetic acid iso-propyl ester

SOLVAGREEN® ≥99 %, for synthesis

CAS No.: 108-21-4

print date: 01.01.2022

formula: $C_5H_{10}O_2$

density: 0,87

molecular weight: 102,13 g/mol

Type analysis

Assay ≥99.0 %

Refractive index 1.376-1.379

Density (20 °C) 0.871-0.875 g/cm³

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

Dr. R. Niemand

Head of Quality Assurance

S. Lorsee

Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5

76185 Karlsruhe

Telefon 0721/5606-0

Telefax 0721/5606-149

E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

