

Kyselina askorbová

Ascorbic acid

Vitamin C

$C_6H_8O_6$

M_r 176,13

CAS: 50-81-7

EINECS: 200-066-2

18650

p.a. / A.G.

Obsah / Assay.....	min. 99,5%
Síranový popel / Sulphated ash	max. 0,05%
Ztráta sušením / Loss on drying (at 105°C)	max. 0,1%
Cl ⁻	max. 0,005%
SO ₄ ²⁻	max. 0,002%
Fe	max. 0,001%
Těžké kovy / Heavy metals (Pb).....	max. 0,001%
pH (5%, 20 °C)	2,2-2,5
[α] _D ²⁰ (c=10 in H ₂ O)	+21 +/- 0,5°