

Specification

Item number: 0752

Methylcyclohexane

≥99 %, for synthesis

CAS No.: 108-87-2

print date: 14.08.2018

formula:

density: 0,77

molecular weight:

Type analysis

Assay (GC)	≥99,0 %
Density (20 °C)	0,768-0,770 g/cm ³
Water	≤0,1 %
Refractive index n_D^{20}	1,420-1,427

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

U. Jeckel
Quality Assurance

S.Hildebrandt
Quality Management

This document is computer printout and has therefore not been signed by hand.