

Specification

Item number:	8131		
Bromide bromate solution			
0,05 mol Br ₂ /l - 0,1 N volumetric standard solution			
CAS No.:		print date:	02.02.2018
formula:			
density:			
molecular weight:			

Type analysis

Amount-of-substance concentration (20 °C)	0,05 mol Br ₂ /l ±0,2 %
Titer	0,998-1,002
Composition:	Potassium bromate, Potassium bromide and Water
Indicator:	Starch
Reference material:	indirectly Potassium dichromate

Reference material is traceable to NIST.

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

U. Jeckel
Quality Assurance

S.Hildebrandt
Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Vorsitzender des Aufsichtsrats: Eberhard Gaul
Geschäftsführer: Lothar Haidmann

