

Specification

Item number: 2957

ROTI®DipSlide CLED/MacCon

ISO 11133, ready-to-use, sterile, for microbiology

CAS No.:

print date:

05.09.2024

formula:

density:

storage temperature: +15 bis +25 °C

molecular weight:

Composition

Agar side 1:

Casein peptone	4 g/l
Gelatine peptone	4 g/l
Beef extract	3 g/l
Lactose	10 g/l
L-Cystine	0.128 g/l
Bromothymol blue	0.02 g/l
Agar	15 g/l
pH-value	7.3±0.2
Colour	green
Growth of test strains	complies

Agar side 2:

Gelatin peptone (Pancreas hydrolysate)	17 g/l
Meat and casein peptone	3 g/l
Lactose monohydrate	10 g/l
Bile salt (Cholate)	1.5 g/l
Sodium chloride (NaCl)	5 g/l
Neutral red	0.03 g/l
Crystal Violet	0.001 g/l
Agar	13.5 g/l
pH-value	7.1±0.2
Colour	violet-red
Growth of test strains	complies

Specification



Detection of bact. of the urinary tract and Enterob./Coliforms

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

Dr. R. Niemand
Head of Quality Assurance

S. Lorse
Head of Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Geschäftsführer: André Houdelet

