



Specification

Item number: 2848

Roti®-Mount Aqua

for microscopy

CAS No.:

print date:

30.05.2018

formula:

density: 1,05-1,18

molecular weight:

Type analysis

Appearance

colourless, viscous liquid

Viscosity (20 °C)

170-220 mPas

Refractive index (20 °C)

1,3750-1,3950

Density (20 °C)

1,05-1,18 g/cm3

Our products are tested for laboratory use only.

The stated information refers to the current status of product quality.

We reserve the right to implement necessary changes.

U. Jeckel

Quality Assurance

S.Hildebrandt

Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Vorsitzender des Aufsichtsrats: Eberhard Gaul
Geschäftsführer: Lothar Haidmann

