

Specification

Item number:	CP62		
Bile Esculin Azide Agar			
DEV, ISO 7899-2:2000, for microbiology			
CAS No.:		print date:	20.04.2017
formula:			
density:		storage temperature:	+20 °C
molecular weight:			

Guarantee analysis

Tryptone	17 g/l
Peptone	3 g/l
Yeast extract	5 g/l
Ox bile dehydrated	10 g/l
Sodium chloride (NaCl)	5 g/l
Esculin	1 g/l
Ammonium iron(III)-citrate	0,5 g/l
Sodium azide	0,15 g/l
Agar	15 g/l
pH-value	7,1 ±0,1

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

U. Jeckel
Quality Assurance

S.Hildebrandt
Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Vorsitzender des Aufsichtsrats: Eberhard Gaul
Geschäftsführer: Lothar Haidmann

