

VELP Scientifica Solutions for Stirring



Stirring speed rpm	up to 1100	up to 1500	up to 1100	up to 1100	up to 1200	up to 1200	80 to 1500	80 to 1500	up to 1100	up to 1100	
Stirring volume L	up to 5	up to 5	up to 5	up to 5	up to 8	up to 25	6 position: up to 2.4 (400 ml each) 15 position: up to 3.75 (250 ml each)	6 position: up to 2.4 (400 ml each) 15 position: up to 3.75 (250 ml each)	up to 5	up to 20	
Timer		•						•			
Reverse Operation		•						•			
Dimensions (WxHxD) mm	120x50x145	120x50x145	120x48x128	160x33x230	171x75x190	250x120x285	230x51.5x370	230x51.5x370	150x55x270	600x55x270	
(in)	(4,7x1,97x5,7)	(4,7x1,97x5,7)	(4.7x1.9x5.0)	(6.3x1.3x9.0)	(6.7x2.9x7.5)	(9.8x4.7x11.2)	(9.0x2.0x14.5)	(9.0x2.0x14.5)	(5.9x2.2x10.6)	(23.6x2.2x10.6)	
Weight	Kg	0.4	0,6	0.55	0.9	1.8	3.7	1.75	2.1	1.2	4.0
	(lb)	(0.9)	-1,32	(1.1)	(2.0)	(4.0)	(8.1)	(3.8)	(4.6)	(2.6)	(8.8)
Power Supply	100÷240 V	100÷240 V	100÷240 V	100÷240 V	115 or 230 V	115 or 230 V	100÷240 V	100÷240 V	100÷240 V	100÷240 V	
Power	0.6 W	0.6 W	0.6 W	5 W	40 W	15 W	3.6 W	9 W	1.2 W	4.8 W	