

Specification

Item number: HN00

Sodium chloride

≥99.8 %, CELLPURE®

CAS No.: 7647-14-5

print date: 19.04.2017

formula: NaCl

density: 2,17

molecular weight: 58,44 g/mol

Guarantee analysis

Endotoxins	≤5 I.U./g
Assay (dry substance)	≥99,8 %
Appearance of solution	complies
Loss on drying	≤0,1 %
Bromide (Br)	≤0,005 %
Iodide (I)	complies
Sulphate (SO ₄)	≤0,02%
Phosphate (PO ₄)	≤0,002 %
Calcium, magnesium, alkaline-earth metals (as Ca)	≤0,01 %
Heavy metals (as Pb)	≤0,0005 %
Aluminium (Al)	≤0,00002 %
Arsenic (As)	≤0,0001 %
Iron (Fe)	≤0,0002 %

Our products are tested for laboratory use only.
The stated information refers to the current status of product quality.
We reserve the right to implement necessary changes.

U. Jeckel
Quality Assurance

S.Hildebrandt
Quality Management

This document is computer printout and has therefore not been signed by hand.

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe

Telefon 0721/5606-0
Telefax 0721/5606-149
E-Mail: info@carlroth.de

Die Firma ist eine Kommanditgesellschaft mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRA 100055. Persönlich haftende Gesellschafterin ist die Firma Roth Chemie GmbH mit Sitz in Karlsruhe, Reg. Gericht Mannheim HRB 100428. Vorsitzender des Aufsichtsrats: Eberhard Gaul
Geschäftsführer: Lothar Haidmann

