

GLYCYL-GLYCIN PUFFERAN®

diglycin, Gly-gly


CAS: 556-50-3

Vzorec: C₄H₈N₂O₃

Technické parametry

Molekulová hmotnost:

132,12



Čistota	Balení	Katalogové číslo
min 99 %, ultra quality	10 g	G 08501
min 99 %, ultra quality	100 g	G 08502
min 99 %, ultra quality	250 g	G 08503
min 99 %, ultra quality	500 g	G 08504
min 98,5 %, Buffer Grade	50 g	R.HN74.1
min 98,5 %, Buffer Grade	100 g	R.HN74.2
min 98,5 %, Buffer Grade	250 g	R.HN74.3
min 98,5 %, Buffer Grade	500 g	R.HN74.4