

Standardy indexu lomu

Standardy pro refraktometrii ROTI®Calipure

Technické parametry

Přesnost: $\pm 0,005$ jednotek

Obal: sklo



Index lomu (20 °C)	Balení	Katalogové číslo
1,38779	15 ml	R.3100.1
1,40485	15 ml	R.3102.1
1,42345	15 ml	R.3103.1
1,44468	15 ml	R.3106.1
1,46768	15 ml	R.3120.1
1,49672	15 ml	R.3124.1
1,50044	15 ml	R.3125.1
1,51726	15 ml	R.3126.1
1,53660	15 ml	R.3143.1
1,65808	15 ml	R.3152.1