

TRIS HYDROCHLORID PUFFERAN®

hydrochlorid tris-(hydroxymethyl)-aminomethanu

CAS: 1185-53-1

Vzorec: $C_4H_{11}NO_3 \cdot HCl$


Technické parametry

Molekulová hmotnost:

157,6

Bezp. věty (GHS):

H315-H319-H335



Čistota	Balení	Katalogové číslo
min 99 %, p.a., PUFFERAN®	50 g	R 90901
min 99 %, p.a., PUFFERAN®	250 g	R.9090.6
min 99 %, p.a., PUFFERAN®	500 g	R 90902
min 99 %, p.a., PUFFERAN®	1 kg	R 90903
min 99 %, p.a., PUFFERAN®	2,5 kg	R.9090.5
min 99 %, p.a., PUFFERAN®	5 kg	R 90904